


Phothera

400 LT

The Phothera 400 LT is a compact full-body phototherapy device designed to treat widespread conditions in space-constrained environments.

At just four feet tall, its lightweight design is easy to maneuver and provides broad treatment coverage in a smaller/shorter footprint. The 400 LT is integrated with the ClearLink® Controller with Guided Mode technology to support a safe and convenient home treatment.

1. The Phothera 400 LT is made with 4 lamps delivering NB-UVB.
2. The Phothera 400 LT features a slim, lightweight, 4-ft design that is easy to set up and maneuver, while its wall-mountable configuration helps maximize space in smaller environments such as dorms, trailers, and apartments.
3. With (4) four foot lamps, advanced ballasts, and precision-formed reflectors, the Phothera 400 LT delivers higher output for shorter treatment times.
4. The Phothera 400 LT offers simplicity and ease of use, featuring an integrated touchscreen controller with Guided Mode technology to support patient-directed home phototherapy and consistent treatment.



TREATMENT AREA	19" × 48" 48.3 cm × 121.9 cm
WIDTH (INCL. FEET)	23.5" 59.7 cm
HEIGHT	50.75" 128.9 cm
DEPTH (INCL. FEET)	16" 40.6 cm
WEIGHT	33 lb 15 kg

For more information, visit phothera.com or contact us at info@phothera.com.