



Phothera

600 LT

The Phothera 600 LT is a full-body phototherapy device used to treat widespread conditions.

Featuring a lightweight, space-saving design, it is equipped with four NB-UVB lamps and is integrated with ClearLink® Controller with Guided Mode technology to support safe and convenient home treatment.

- 1. The Phothera 600 LT is made with 4 lamps delivering NB-UVB.
- 2. The Phothera 600 LT features a slim, lightweight design with easy setup, making it quick to assemble, move into place, and begin treatment.
- 3. With (4) six foot lamps, advanced ballasts, and precision-formed reflectors, the Phothera 600 LT delivers higher output for shorter treatment times.
- 4. The Phothera 600 LT offers simplicity and ease of use, featuring an integrated touchscreen controller with Guided Mode technology to support patient-directed home phototherapy and consistent treatment.



TREATMENT AREA	19" x 70" 48.3 cm x 177.8 cm
WIDTH (INCL. FEET)	23.5" 59.7 cm
HEIGHT	72.75" 184.8 cm
DEPTH (INCL. FEET)	16" 40.6 cm
WEIGHT	40 lb 18 kg

For more information, visit phothera.com or contact us at info@phothera.com.